

SMILE – you're on candid camera

Imaging and AI in zooplankton studies

Department of Estuarine & Delta Systems, NIOZ

Pieter Hovenkamp
Dick van Oevelen
Karline Soetaert

Wageningen Marine Research


Lodewijk van Walraven
Jeroen Hoekendijk
Robbert Jak

University College Roosevelt (UU/UCR)


Frank van der Stappen



Zooplankton in the North Sea food web






Multitrophic pelagic North Sea food web



- Important link in the food web
- Biological carbon pump
- Many benthic organisms have a larval stage

Sampling is limited (costly) and methodologically difficult

Continuous plankton recorder (CPR)



- No depth info
- Time-consuming
- No coastal data
- Sampling bias
- Gelatinous taxa missing



Three plankton imaging case studies:

- In situ imaging – Doggersbank
- On-board imaging – Coastal North Sea
- Mobile imaging – in development




In situ imaging – Doggersbank





8.83 cm



Deployment time: 0:00:00

Number of frames: 1

Megabytes: 24.5

Number of image segments: 90








Figure from Kröncke & Knust, 1995


72 million segments



Sharp chlorophyll maximum



Mixed chlorophyll-a layer




Three plankton imaging case studies:

- In situ imaging – Doggersbank
- On-board imaging – Coastal North Sea
- Mobile imaging – in development



On-board imaging – Coastal North Sea


- Plankton Imager (PI10) developed by Plankton Analytics with CEFAS
- Autonomous and continuous underway sampling of mesoplankton composition
- Does not interfere with ship operations



On-board imaging – Coastal North Sea




- MONS cruise with R/V Tridens
- Continuous sampling along the transect
 - 149 hours
 - 980 nm
 - 1.125 billion particles detected
 - 86 million images stored (7.6 %)



RV Zirfaea
(water quality)




RV Tridens
(fisheries)




On-board imaging – Coastal North Sea


Holoplankton (incl. gelatinous)




Meroplankton




Copepods




Oikopleura




Phaeocystis



Exuvia (empty skins)



Noctiluca



All Larvaceans: >40 ind L⁻¹



Three plankton imaging case studies:

- In situ imaging – Doggersbank
- On-board imaging – Coastal North Sea
- Mobile imaging – in development



Mobile imaging – in development


REDstack:


Company that develops membrane stacks to generate electricity by Reverse ElectroDialysis (RED)

For upscaling, large (up to 100m³/s) amounts of seawater are needed

Requirements:

- Mortality of plankton (i.e. copepods and (fish) larvae) needs to be reduced to avoid ecological impact
- Collection of water with fine particles (>10µm) needs to be reduced to avoid pollution






Development of a mobile container for mobile and real-time monitoring

Check out the poster of Amin Niamir

Measurement of zooplankton and particles at facility and in Wadden sea



Summary

Copepods

